

Corus CAD is the only clinically validated blood test that helps clinicians in the office setting evaluate whether a patient's symptoms are due to obstructive coronary artery disease (CAD).^{*} Genes that are expressed in cells of the innate and adaptive immune systems were examined during the development of Corus CAD. The algorithm genes reflect cellular and molecular processes associated with:


- Cellular apoptosis and necrosis
- Inflammatory cell migration into atherosclerotic plaque
- Innate and adaptive immune response to LDL oxidation and other inflammatory processes

The Corus CAD algorithm incorporates age, sex, and gene expression and was validated in the PREDICT and COMPASS multi-center studies.^{1,2}

Algorithm Terms

The genes in the Corus CAD algorithm are grouped into six terms based on correlation of gene expression. The Corus CAD algorithm integrates age, sex, and gene expression to calculate a score. The complete algorithm generates a score that helps inform the clinician on the likelihood of obstructive CAD.

FIGURE 1: CORUS CAD GENE TERM SCHEMATIC



Term 1

This term contains genes that likely reflect increased neutrophil apoptosis (*CASP5*, *TNFRSF10C*),^{3,4} activation of the interleukin-18 and anti-inflammatory signaling pathways (*IL18RAP*, *TNFAIP6*),^{5,6} and other genes involved in the innate immune response and the recruitment of neutrophils to sites of inflammation (*TLR4*, *IL8RB/CXCR2*).^{7,8,9}

Term 2

The main component of this term is comprised of three genes primarily found in granulocytes and neutrophils. Two of the genes (*S100A12*, *S100A8*)^{10,11} have been implicated in coronary disease and cardiovascular events, and the third (*CLEC4E*)¹² has been shown to play a role in the surveillance of cellular necrosis.

[Continued on reverse]

^{*}Obstructive coronary artery disease (CAD) is defined as at least one atherosclerotic plaque causing $\geq 50\%$ luminal diameter stenosis in a major coronary artery ($\geq 1.5\text{mm}$ lumen diameter) as determined by invasive quantitative coronary angiography (QCA) or core-lab computerized tomography angiography (CTA) ($\geq 2.0\text{mm}$ lumen diameter).

Term 3

This term is comprised of two pairs of genes reflecting natural killer cell activation (SLAMF7, KLRC4)^{13,14} normalized to genes expressed in T lymphocytes (TMC8, CD3D).^{15,16}

Term 4

This term consists of two groups of genes reflecting levels of B (SPIB, CD79B)^{17,18} and T lymphocytes (TMC8, CD3D). The B to T cell ratio is likely sensitive to the presence of atherosclerosis.

Term 5

This term includes AF289562, a gene for which the underlying biology has yet to be determined. AF289562 is normalized to two genes (TFCP2, HNRPF) that are cell-type independent in their blood expression.

Term 6

This term consists of TSPAN16,¹⁵ a gene for which the exact mechanism and associated biological pathways are unknown. TSPAN16 is normalized to two genes (TFCP2, HNRPF) that are cell-type independent in their blood expression. This term is used only in males.

For the Corus CAD Intended Use Statement please refer to www.cardiodx.com/corus-cad-intended-use

REFERENCES:

1. Rosenberg S, Elashoff MR, Beineke P, et al. Multicenter Validation of the Diagnostic Accuracy of a Blood-Based Gene Expression Test for Assessing Obstructive Coronary Artery Disease in Nondiabetic Patients. *Ann Intern Med.* 2010;153:425-434.
2. Thomas GS, Voros S, McPherson JA, et al. A Blood-Based Gene Expression Test for Obstructive Coronary Artery Disease Tested in Symptomatic Nondiabetic Patients Referred for Myocardial Perfusion Imaging: The COMPASS Study. *Circ Cardiovasc Genet.* 2013;6(2):154-162.
3. Pop C, Salvesen GS. Human Caspases: Activation, Specificity, and Regulation. *J Biol Chem.* 2009;284(33):21777-81.
4. Sheridan JP, Marsters SA, Pitti RM, et al. Control of TRAIL-Induced Apoptosis by a Family of Signaling and Decoy Receptors. *Science.* 1997;277(5327):818-21.
5. Tired L, Godefroy T, Lubos E, et al. Genetic Analysis of the Interleukin-18 System Highlights the Role of the Interleukin-18 Gene in Cardiovascular Disease. *Circulation.* 2005;112(5):643-50.
6. Milner CM, Higman VA, Day AJ. TSG-6: A Pluripotent Inflammatory Mediator? *Biochem Soc Trans.* 2006;34(Pt 3):446-50.
7. de Kleijn D, Pasterkamp G. Toll-Like Receptors in Cardiovascular Diseases. *Cardiovasc Res.* 2003;60(1):58-67.
8. den Dekker WK, Cheng C, Pasterkamp G, et al. Toll-Like Receptor 4 in Atherosclerosis and Plaque Destabilization. *Atherosclerosis.* 2010;209(2):314-20.
9. Bernhagen J, Krohn R, Lue H, et al. MIF is a Noncognate Ligand of CXC Chemokine Receptors in Inflammatory and Atherogenic Cell Recruitment. *Nat Med.* 2007;13(5):587-96.
10. Mahajan N, Malik N, Bahl A, et al. Receptor for Advanced Glycation End Products (RAGE) and its Inflammatory Ligand EN-RAGE in Non-Diabetic Subjects with Pre-Mature Coronary Artery Disease. *Atherosclerosis.* 2009;207(2):597-602.
11. Ionita MG, Vink A, Dijke IE, et al. High Levels of Myeloid-Related Protein 14 in Human Atherosclerotic Plaques Correlate With the Characteristics of Rupture-Prone Lesions. *Arterioscler Thromb Vasc Biol.* 2009;29(8):1220-7.
12. Yamasaki S, Ishikawa E, Sakuma M, et al. Mincle is an ITAM-Coupled Activating Receptor That Senses Damaged Cells. *Nat Immunol.* 2008;9(10):1179-1188.
13. Cruz-Munoz ME, Dong Z, Shi X, et al. Influence of CRACC, a SLAM Family Receptor Coupled to the Adaptor EAT-2, on Natural Killer Cell Function. *Nat Immunol.* 2009;10(3):297-305.
14. Kim DK, Kabat J, Borrego F, et al. Human NKG2F is Expressed and can Associate with DAP12. *Mol Immunol.* 2004;41(1):53-62.
15. Wu C, Orozco C, Boyer J, et al. BioGPS: An Extensible and Customizable Portal for Querying and Organizing Gene Annotation Resources. *Genome Biol.* 2009;10(11):R130.
16. Wang M, Windgassen D, Papoutsakis ET. Comparative Analysis of Transcriptional Profiling of CD3+, CD4+ and CD8+ T Cells Identifies Novel Immune Response Players in T-Cell Activation. *BMC Genomics.* 2008;9:225.
17. Garrett-Sinha LA, Su GH, Rao S, et al. PU.1 and Spi-B are Required for Normal B Cell Receptor-Mediated Signal Transduction. *Immunity.* 1999;10(4):399-408.
18. Vasile S, Coligan JE, Yoshida M, et al. Isolation and Chemical Characterization of the Human B29 and MB-1 Proteins of the B Cell Antigen Receptor Complex. *Mol Immunol.* 1994;31(6):419-27.

CardioDx

600 Saginaw Drive
Redwood City, CA 94063

www.cardiodx.com

Customer Service

service@cardiodx.com
P 866.941.4996

Medical Affairs

medicalaffairs@cardiodx.com
P 866.941.4996

